

A woman with short blonde hair, wearing a light green t-shirt and white shorts, stands on a stone path in a garden. She is holding a green watering can in her left hand and a pink flower in her right hand. She is wearing a grey knee brace on her right leg. The garden is filled with various plants, including tall grasses and pink flowers. The background shows a clear blue sky and more greenery.

**THERE'S
ONLY ONE.**

THE UNLOADER ONE® BY ÖSSUR



Life Without Limitations®

UNLOADER ONE.[®] THE ORIGINAL. THE STANDARD. THE UNLOADER.

3 cheers for 3 points.

The efficacy of our 3-Point Leverage System has been clinically proven by multiple peer-reviewed studies. Unloader One takes our patented GII[®] technology to the next level.



Ratchet it up... or down.

Featuring a tool-free ratchet, Unloader One enables on-the-fly adjustable and repeatable unloading, helping patients better manage their unicompartmental OA knee pain.

Barely there.

As the lightest and lowest-profile OA knee brace on the market, Unloader One is barely noticeable to the user and fits inconspicuously and comfortably under clothing.



Migration is for the birds.

Other OA knee braces often migrate. Össur's patented Sensil[®] silicone liners are specially designed to comfortably secure the Unloader One in place, maximizing effectiveness.

30 DAY
GUARANTEE

Prescribe the Unloader One for any patient. If they are not happy with it for whatever reason, they can return it within 30 days for a full refund.



CLINICALLY PROVEN OSTEOARTHRITIS PAIN RELIEF



RECENT CLINICAL RESEARCH BY THE STEADMAN PHILIPPON RESEARCH INSTITUTE* ESTABLISHED THAT THE UNLOADER ONE® DECREASES PAIN, IMPROVES FUNCTION AND REDUCES THE USAGE OF PAIN MEDICATION.

STEADMAN PHILIPPON RESEARCH RESULTS

Patients reported a significant improvement in disability (per WOMAC score):

Knee compartment unloaded	WOMAC Score improvement
Medial	16 points
Lateral	8 points

At 6 months, a significant percentage of patients reported a reduction in medication use:

% of patients reporting medication reduction	Type of medication
27%	Pain medication
31%	OTC anti-inflammatories
35%	Prescription anti-inflammatories

Median patient satisfaction: 8 (out of 10)

* Research available at www.ossur.com/unloader

EASY TO FIT AND WEAR



Alligator clips make strapping reversible for external or internal fastening



Pre-attached popliteal pad improves patient comfort and maintains proper strap positioning



Color-coded buckles for easy and intuitive brace application



Matching dots on straps enable easier re-assembly of brace after washing



Life Without Limitations®

AVAILABLE IN PLUS SIZES, JUST LIKE YOUR PATIENTS



Comfortable containment.

For a variety of reasons, many sufferers of OA knee pain are overweight. For these patients, containing excess soft tissue is essential if they are to be braced effectively and comfortably. That's why our new Unloader One Plus model is equipped with a wider thigh strap¹ than the standard version, plus a gastroc strap that helps secure the brace in place.²

Extra large support.

Bigger legs require bigger support. So, for the new Unloader One Plus model, we included a sturdier upright with an Adjustable Dynamic Joint that provides clinicians with additional leverage for optimal pain management. Multiple, peer-reviewed clinical and biomechanical studies have proven the efficacy of this technology.



UNLOADER ONE® ORDERING INFORMATION

UNLOADER ONE® CUSTOM

Color	Right	Left
Grey	B-115501711	B-115601711
Black + 12 Colors (see below)	B-115500111	B-115600111

UNLOADER ONE® OTS

Size	Leg Circumference*	Medial	Lateral
S	12.2-14" (31-35.6cm)	B-240X19712	B-240X29712
M	14-16.75" (35.6-42.5cm)	B-240X19713	B-240X29713
L	16.75-20" (42.5-50.8cm)	B-240X19714	B-240X29714
XL	20-24" (50.8-61cm)	B-240X19715	B-240X29715
* 6" below mid-patella		X = 5 (Right), 6 (Left)	

Short Version

Size	Leg Circumference*	Medial	Lateral
XS	9.8-12.2" (24.9-31cm)	B-241X19711	B-241X29711
S	12.2-14" (31-35.6cm)	B-241X19712	B-241X29712
M	14-16.75" (35.6-42.5cm)	B-241X19713	B-241X29713
* 6" below mid-patella		X = 5 (Right), 6 (Left)	

Notes:

Accessories available: Straps, pads, liners, AMS wrap, buckles, alligator tabs
Short Length: 14.5" (36.8cm); Standard Length: 15.5" (39.4cm); Long Length: 16.5" (41.9cm)
No weight limit. Sizing based on thigh/calf circumference.

UNLOADER ONE® PLUS CUSTOM

Color	Right	Left
Grey	B-119501700	B-119601700
Black + 12 Colors (see below)	B-119500100	B-119600100

UNLOADER ONE® PLUS OTS

Size	Thigh Circ*	Calf Circ*	Part Number
L	21-25" (53-64cm)	14-16.75" (36-43cm)	B-245XY0004
XL	25-29" (64-74cm)	16.75-20" (43-51cm)	B-245XY0005
* 6" above/below mid-patella			X = 5 (Right) or 6 (Left) Y = 1 (Medial) or 2 (Lateral)

Long Version

Size	Thigh Circ*	Calf Circ*	Part Number
L	21-25" (53-64cm)	14-16.75" (36-43cm)	B-244XY0004
XL	25-29" (64-74cm)	16.75-20" (43-51cm)	B-244XY0005
* 6" above/below mid-patella			X = 5 (Right) or 6 (Left) Y = 1 (Medial) or 2 (Lateral)

Custom Colors

SB	OB	NB	GR	CH	RD	PK	WH	SV	CC	YW	BK
----	----	----	----	----	----	----	----	----	----	----	----

SB = Sky Blue, OB = Ocean Blue, NB = Navy Blue, GR = Green, CH = Champagne, RD = Red, PK = Pink, WH = White, SV = Silver, CC = Charcoal, YW = Yellow, BK = Black

USA 800 233 6263
CANADA 800 663 5982
WEB WWW.OSSUR.COM

FOLLOW ÖSSUR ON

